

Product Datasheet

Mouse CHD4 ELISA Kit (orb438613)

Description	Mouse Chromodomain Helicase DNA Binding Protein 4 (CHD4) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q14839
Sensitivity	0.058ng/mL
Entrez	1108
Expiration Date	Please enquire.