

Product Datasheet

Human CHGA ELISA Kit (orb438618)

Description	Human Chromogranin A (CHGA) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P10645
Sensitivity	23.5pg/mL
Entrez	1113
Expiration Date	Please enquire.