

Product Datasheet

Mouse CHRM2 ELISA Kit (orb438630)

Description	Mouse Cholinergic Receptor, Muscarinic 2 (CHRM2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9ERZ4
Sensitivity	0.057ng/mL
Entrez	243764
Expiration Date	Please enquire.