

Product Datasheet

Rat CHRM2 ELISA Kit (orb438631)

Description	Rat Cholinergic Receptor, Muscarinic 2 (CHRM2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P10980
Sensitivity	0.057ng/mL
Entrez	81645
Expiration Date	Please enquire.