

Product Datasheet

Mouse CKB ELISA Kit (orb438647)

Description	Mouse Creatine Kinase, Brain (CKB) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q04447
Sensitivity	6.9pg/mL
Entrez	12709
Expiration Date	Please enquire.