

Product Datasheet

Bovine CKM ELISA Kit (orb438648)

Description	Bovine Creatine Kinase, Muscle (CKM) Sandwich ELISA Kit
Reactivity	Bovine
Range	12.5-800U/L
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9XSC6
Sensitivity	4.65U/L
Entrez	286822
Expiration Date	Please enquire.