

Product Datasheet

Canine CKM ELISA Kit (orb438652)

Description	Canine Creatine Kinase, Muscle (CKM) Sandwich ELISA Kit
Reactivity	Canine
Range	6.25-400U/L
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P05123
Sensitivity	3.1U/L
Entrez	476435
Expiration Date	Please enquire.