

Product Datasheet

Human CKMT1A ELISA Kit (orb438659)

Description	Human Creatine Kinase, Mitochondrial 1A (CKMT1A) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P12532
Sensitivity	0.57ng/mL
Entrez	548596
Expiration Date	Please enquire.