

Product Datasheet

Rat CLDN3 ELISA Kit (orb438672)

Description	Rat Claudin 3 (CLDN3) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q63400
Sensitivity	11.5pg/mL
Entrez	65130
Expiration Date	Please enquire.