

Product Datasheet

Human CLEC11A ELISA Kit (orb438680)

Description	Human C-Type Lectin Domain Family 11, Member A (CLEC11A) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9Y240
Sensitivity	0.51ng/mL
Entrez	6320
Expiration Date	Please enquire.