

Product Datasheet

Human CLEC3B ELISA Kit (orb438682)

Description	Human C-Type Lectin Domain Family 3, Member B (CLEC3B) Sandwich ELISA Kit
Reactivity	Human
Range	0.187-12ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	P05452
Sensitivity	0.07ng/mL
Entrez	7123
Expiration Date	Please enquire.