

Product Datasheet

Human CLIC4 ELISA Kit (orb438685)

Description	Human Chloride Intracellular Channel Protein 4 (CLIC4) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9Y696
Sensitivity	30pg/mL
Entrez	25932
Expiration Date	Please enquire.