

Product Datasheet

Human CLMP ELISA Kit (orb438687)

Description	Human Coxsackie And Adenovirus Receptor Like Membrane Protein (CLMP) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9H6B4
Sensitivity	0.058ng/mL
Entrez	79827
Expiration Date	Please enquire.