

Product Datasheet

Human CMA1 ELISA Kit (orb438694)

Description	Human Chymase 1, Mast Cell (CMA1) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P23946
Sensitivity	6.5pg/mL
Entrez	1215
Expiration Date	Please enquire.