

Product Datasheet

Mouse CMA1 ELISA Kit (orb438695)

Description	Mouse Chymase 1, Mast Cell (CMA1) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P21844
Sensitivity	30.15pg/mL
Entrez	17228
Expiration Date	Please enquire.