

Product Datasheet

Rat CMA1 ELISA Kit (orb438696)

Description	Rat Chymase 1, Mast Cell (CMA1) Sandwich ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P50339
Sensitivity	0.57ng/mL
Entrez	25627
Expiration Date	Please enquire.