

Product Datasheet

Human CNP ELISA Kit (orb438700)

Description	Human C-Type Natriuretic Peptide (CNP) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P23582
Sensitivity	26.4pg/mL
Entrez	4880
Expiration Date	Please enquire.