

Product Datasheet

Rat CNP ELISA Kit (orb438702)

Description	Rat C-Type Natriuretic Peptide (CNP) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P13233
Sensitivity	12.4pg/mL
Entrez	25275
Expiration Date	Please enquire.