

Product Datasheet

Chicken CNTF ELISA Kit (orb438707)

Description	Chicken Ciliary Neurotrophic Factor (CNTF) Sandwich ELISA Kit
Reactivity	Gallus
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q02011
Sensitivity	29.4pg/mL
Entrez	396488
Expiration Date	Please enquire.