

Product Datasheet

Porcine CNTF ELISA Kit (orb438710)

Description	Porcine Ciliary Neurotrophic Factor (CNTF) Sandwich ELISA Kit
Reactivity	Porcine
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	F1RMF2
Sensitivity	22.4pg/mL
Entrez	100518821
Expiration Date	Please enquire.