

Product Datasheet

Rat CNTF ELISA Kit (orb438711)

Description	Rat Ciliary Neurotrophic Factor (CNTF) Sandwich ELISA Kit
Reactivity	Rat
Range	39.062-2500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P20294
Sensitivity	17.5pg/mL
Entrez	25707
Expiration Date	Please enquire.