

Product Datasheet

Rabbit CNTF ELISA Kit (orb438712)

Description	Rabbit Ciliary Neurotrophic Factor (CNTF) Sandwich ELISA Kit
Reactivity	Rabbit
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P14188
Sensitivity	11.8pg/mL
Entrez	100009117
Expiration Date	Please enquire.