

Product Datasheet

Human COL11a1 ELISA Kit (orb438720)

Description	Human Collagen Type XI Alpha 1 (COL11a1) Sandwich ELISA Kit
Reactivity	Human
Range	39.062-2500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P12107
Sensitivity	15pg/mL
Entrez	1301
Expiration Date	Please enquire.