

Product Datasheet

Rat COL19 ELISA Kit (orb438729)

Description	Rat Collagen Type XIX (COL19) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	D3ZCQ0
Sensitivity	0.116ng/mL
Entrez	367236
Expiration Date	Please enquire.