


Product Datasheet

H-Ras antibody (orb43873)

Description	Mouse Monoclonal to H-Ras.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to amino acids DIHQYREQIKRVKDSDDC of human H-Ras protein
Target	H-RAS
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Western blotting: Recommended dilution: 1-2 µg/ml, reducing conditions.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P01112
Entrez	3265
Dilution Range	Western blotting: Recommended dilution: 1-2 µg/ml, reducing conditions.
Expiration Date	12 months from date of receipt.



Western blotting analysis of human H-Ras...