

## Product Datasheet

### Human COL1a1 ELISA Kit (orb438730)

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Collagen Type I Alpha 1 (COL1a1)<br>Sandwich ELISA Kit                             |
| <b>Reactivity</b>        | Human                                                                                    |
| <b>Range</b>             | 0.312-20ng/mL                                                                            |
| <b>Form/Appearance</b>   | Traditional                                                                              |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                    |
| <b>Application notes</b> | serum, plasma and other biological fluids.                                               |
| <b>Sample Types</b>      | serum, plasma and other biological fluids.                                               |
| <b>Uniprot ID</b>        | <b>P02452</b>                                                                            |
| <b>Sensitivity</b>       | 0.133ng/mL                                                                               |
| <b>Entrez</b>            | <b>1277</b>                                                                              |
| <b>Expiration Date</b>   | Please enquire.                                                                          |