

Product Datasheet

Mouse COL1a1 ELISA Kit (orb438731)

Description	Mouse Collagen Type I Alpha 1 (COL1a1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P11087
Sensitivity	0.137ng/mL
Entrez	12842
Expiration Date	Please enquire.