

Product Datasheet

Rabbit COL1a2 ELISA Kit (orb438738)

Description	Rabbit Collagen Type I Alpha 2 (COL1a2) Sandwich ELISA Kit
Reactivity	Rabbit
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q28668
Sensitivity	0.104ng/mL
Entrez	100008997
Expiration Date	Please enquire.