

Product Datasheet

Human COL4a1 ELISA Kit (orb438754)

Description	Human Collagen Type IV Alpha 1 (COL4a1) Sandwich ELISA Kit
Reactivity	Human
Range	0.937-60ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P02462
Sensitivity	0.45ng/mL
Entrez	1282
Expiration Date	Please enquire.