

Product Datasheet

Mouse COL4a1 ELISA Kit (orb438755)

Description	Mouse Collagen Type IV Alpha 1 (COL4a1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P02463
Sensitivity	0.056ng/mL
Entrez	12826
Expiration Date	Please enquire.