

Product Datasheet

Mouse COL5a2 ELISA Kit (orb438762)

Description	Mouse Collagen Type V Alpha 2 (COL5a2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q3U962
Sensitivity	0.116ng/mL
Entrez	12832
Expiration Date	Please enquire.