

Product Datasheet

Human COL6a1 ELISA Kit (orb438763)

Description	Human Collagen Type VI Alpha 1 (COL6a1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P12109
Sensitivity	0.063ng/mL
Entrez	1291
Expiration Date	Please enquire.