

Product Datasheet

Mouse COL6a1 ELISA Kit (orb438764)

Description	Mouse Collagen Type VI Alpha 1 (COL6a1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q04857
Sensitivity	0.268ng/mL
Entrez	12833
Expiration Date	Please enquire.