

Product Datasheet

Rat COL6a1 ELISA Kit (orb438765)

Description	Rat Collagen Type VI Alpha 1 (COL6a1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	D3ZUL3
Sensitivity	0.062ng/mL
Entrez	294337
Expiration Date	Please enquire.