

Product Datasheet

Rat COL9a1 ELISA Kit (orb438772)

Description	Rat Collagen Type IX Alpha 1 (COL9a1) Sandwich ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	F1LQ93
Sensitivity	0.52ng/mL
Entrez	305104
Expiration Date	Please enquire.