

Product Datasheet

Transferrin antibody (orb43878)

Description

Mouse Monoclonal to Transferrin.

Reactivity

Human, Porcine, Rabbit

Conjugation

Unconjugated

Tested Applications

ELISA, FA, ICC, IHC-P, RIA, WB

Immunogen

Purified porcine transferrin.

Target

Transferrin

Preservatives

Phosphate buffered saline (PBS), pH 7.4

Concentration

1 mg/ml

Storage

Store at 2-8°C. Do not freeze.

Note

For research use only

Application notes

Functional application: The antibody HTF-14 blocks binding of transferrin to its receptor.
Immunohistochemistry (paraffin sections):
Recommended dilution: 10 µg/ml; positive tissue: placenta. Western blotting: non-reducing conditions, recommended dilution: 1-2 µg/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity

Purified by protein-A affinity chromatography.

Uniprot ID

P02787

Entrez


7018

Dilution Range


Functional application: The antibody HTF-14 blocks binding of transferrin to its receptor.
Immunohistochemistry (paraffin sections):
Recommended dilution: 10 µg/ml; positive tissue: placenta. Western blotting: non-reducing conditions, recommended dilution: 1-2 µg/ml.

Expiration Date

12 months from date of receipt.



Human transferrin detected by the mouse ...



Immunohistochemistry staining of human p...