

Product Datasheet

Human COPR5 ELISA Kit (orb438781)

Description	Human Cooperator Of PRMT5 (COPR5) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9NQ92
Sensitivity	0.055ng/mL
Entrez	55352
Expiration Date	Please enquire.