

Product Datasheet

Human COPT1 ELISA Kit (orb438782)

Description	Human Copper Transporter 1 (COPT1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O15431
Sensitivity	0.061ng/mL
Entrez	1317
Expiration Date	Please enquire.