

Product Datasheet

General Cort ELISA Kit (orb438786)

Description	General Corticosterone (Cort) competitive inhibition ELISA Kit
Reactivity	All
Range	7.812-500ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, urine and other biological fluids.
Sample Types	serum, plasma, urine and other biological fluids.
Sensitivity	2.62ng/mL
Expiration Date	Please enquire.