



Product Datasheet

Transferrin antibody (orb43879)

Description	Mouse Monoclonal to Transferrin.
Reactivity	Human, Porcine, Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC-P, RIA, WB
Immunogen	Purified porcine transferrin.
Target	Transferrin
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: placenta. Western blotting: non-reducing conditions, recommended dilution: 1-2 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P02787
Entrez	7018
Dilution Range	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: placenta. Western blotting: non-reducing conditions, recommended dilution: 1-2 µg/ml.
Expiration Date	12 months from date of receipt.



Western blotting detection of human tran...



Immunohistochemistry staining of human p...