

Product Datasheet

Rat CPA1 ELISA Kit (orb438794)

Description	Rat Carboxypeptidase A1, Pancreatic (CPA1) Sandwich ELISA Kit
Reactivity	Rat
Range	1.56-100U/L
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P00731
Sensitivity	0.52U/L
Entrez	24269
Expiration Date	Please enquire.