

Product Datasheet

Rat CPA2 ELISA Kit (orb438795)

Description	Rat Carboxypeptidase A2, Pancreatic (CPA2) Sandwich ELISA Kit
Reactivity	Rat
Range	1.56-100U/L
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	G3V976
Sensitivity	0.8U/L
Entrez	296959
Expiration Date	Please enquire.