

Product Datasheet

Human CPA3 ELISA Kit (orb438796)

Description	Human Carboxypeptidase A3, Mast Cell (CPA3) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P15088
Sensitivity	34pg/mL
Entrez	1359
Expiration Date	Please enquire.