

Product Datasheet

Human CPA4 ELISA Kit (orb438799)

Description	Human Carboxypeptidase A4 (CPA4) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9UI42
Sensitivity	0.051ng/mL
Entrez	51200
Expiration Date	Please enquire.