

## Product Datasheet

### Rat C-Peptide ELISA Kit (orb438804)

|                          |                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Rat C-Peptide competitive inhibition ELISA Kit                                                          |
| <b>Reactivity</b>        | Rat                                                                                                     |
| <b>Range</b>             | 125-8000pg/mL                                                                                           |
| <b>Form/Appearance</b>   | Traditional                                                                                             |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt.                |
| <b>Note</b>              | For research use only                                                                                   |
| <b>Application notes</b> | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.   |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| <b>Sensitivity</b>       | 33.27pg/mL                                                                                              |
| <b>Expiration Date</b>   | Please enquire.                                                                                         |