

Product Datasheet

Human CPM ELISA Kit (orb438807)

Description	Human Carboxypeptidase M (CPM) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P14384
Sensitivity	28pg/mL
Entrez	1368
Expiration Date	Please enquire.