

Product Datasheet

Human CPT1B ELISA Kit (orb438822)

Description	Human Carnitine Palmitoyltransferase 1B, Muscle (CPT1B) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q92523
Sensitivity	0.084ng/mL
Entrez	1375
Expiration Date	Please enquire.