

Product Datasheet

Human CPZ ELISA Kit (orb438825)

Description	Human Carboxypeptidase Z (CPZ) Sandwich ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q66K79
Sensitivity	0.226ng/mL
Entrez	8532
Expiration Date	Please enquire.