

Product Datasheet

Rat CR2 ELISA Kit (orb438828)

Description	Rat Complement Receptor 2 (CR2) Sandwich ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	A0A0G2K684
Sensitivity	0.57ng/mL
Entrez	289395
Expiration Date	Please enquire.