

Product Datasheet

Human CRHBP ELISA Kit (orb438839)

Description	Human Corticotropin Releasing Hormone Binding Protein (CRHBP) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P24387
Sensitivity	0.052ng/mL
Entrez	1393
Expiration Date	Please enquire.